We are IntechOpen, the world's leading publisher of Open Access books Built by scientists, for scientists



186,000

200M



Our authors are among the

TOP 1% most cited scientists





WEB OF SCIENCE

Selection of our books indexed in the Book Citation Index in Web of Science™ Core Collection (BKCI)

Interested in publishing with us? Contact book.department@intechopen.com

Numbers displayed above are based on latest data collected. For more information visit www.intechopen.com



The Matrix Metalloproteinases and Cerebral Ischemia

Wan Yang and Guangqin Li Department of Neurology, the First Affiliated Hospital, Chongqing Medical University China

1. Introduction

Matrix metalloproteinases(MMPs) are a family of zinc-containing endopeptidases which can degrade extracellular matrix (ECM) components at physiological PH. Thirty members of MMPs have been found so far ^[1]. They are widely distributed in plants,vertebrates and invertebrates cells. In human body, MMPs are mainly produced by vascular smooth muscle cells, monocytes, endothelial cells and so on. MMPs are synthesized as latent enzymes (zymogens) that are secreted or membrane-associated and must be proteolytically processed to their active form. Additionally, MMPs can be inhibited by endogenous inhibitors (e.g. TIMPs, tissue inhibitors of metalloproteinases).

1.1 Construction and function of matrix metalloproteinases

1.1.1 Construction of MMPs

MMP members have similar structures. They are usually composed of six structural domains with different functions.

1.2 The signal peptide and propeptide domain

Human MMPs(except MMP14) have a signal peptide sequence. The role of the signal peptide is to guide the post-translation substrate to the cytoplasm endoplasmic reticulum. Propeptide domain contains the conserved sequence of Pro-Arg-Cys-Gly-Val/Asn-Pro-Asp(PRCGV/NPD) which is responsible for maintaining the stability of plasminogen. When the propeptide domain is cut off by the exogenous enzymes, MMPs plasminogen can be activated.

1.3 The catalysis domain

There are two zinc ion binding domains and at least one calcium binding domain in the catalysis domain. Of the 2 zinc ion binding domains, one is in the activation center which is responsible for the catalytic process of MMPs. The other one is the structural zinc ion domain. In the catalysis domain, both gelatinase A and B have a insertion sequence of 175 residues. This insertion sequence is the type-fibronectin binding domain. Studies suggest that this domain may be responsible for the integration between gelatinase and its substrate.

1.4 The hinge domain and hemopexin binding domain

The hinge domain is located between the catalysis and the hemopexin binding domian and is linked with the terminal amino acid residues of the hemopexin binding domain by the disulfide bond. The hemopexin domain contains four duplicate sequences which has the weak homology with hemopexin and vitronectin. This domain is considered to be relevant to most MMPs' substrate specificity and also plays an important role in the integration between MMPs and tissue inhibitors of metalloproteinases (TIMPs).

1.5 The transmembrane domain

This domain exists in the carboxyl terminal of the Membrane-type MMPs(MT-MMPs) and can fix MT-MMPs to the cell membrane.

MMP members have different characteristics based on the six structural domains.

2. Function of matrix metalloproteinases

MMPs are the essential enzymes that may play a role in the degradation of ECM in the connective tissue. They can degrade almost all components of ECM and play an important role in various physiological and pathological processes in human body. Under physiological conditions, MMPs participate in the process of tissue demodeling such as wound healing, bone resorption, pregnancy, childbirth and breast atrophy. Normal physiological process depends on the control and coordination between MMPs and TIMPs. When infection or other stimulation occurs, the expression and activation of MMP are out of control, which may lead to the excessive degradation of ECM.

MMPs family can be divided into the following six categories based on the strucuture of their substrate, sequence similarity and characteristics of structural domain.

2.1 Collagenases

MMP-1, MMP-8, MMP-13 and MMP-18 are included. The main feature of these enzymes is to resolve collagen types I, II, III. Collagenase can also resolve some ECM components and other non-extracellular matrix molecules.

2.2 Gelatinases

Gelatinase A(MMP-2) and gelatinase B(MMP-9) are included. They are secreted or membrane-associated and must be proteolytically processed to their active form. MMP-2 can digest gelatin, collagen types IV, V, VIII, X laminin, elastin and fibronectin. The molecular weight of MMP-9 is 92KD which is the largest one. MMP-9 can be synthesized by various cells, such as astrocytes, vascular endothelial cells, microglias, neutrophils and macrophages. MMP-9 is mainly affected by the regulation of plasminogen activator and its subsrates incluce gelatinase, collagen types IV, V and elastin^[2]. Researches have shown that MMP-2 and MMP-9 may play a role in angiogenesis^[3], atherosclerosis^[4,5] and ischemic brain injury^[6].

2.3 Stromelysins

Stromelysin-1(MMP-3) and stromelysin-2(MMP-10) have the same substrate specificity. However, the proteolytic efficiency of MMP-3 is higher than that of MMP-10. In addition to digesting extracellular matrix components, MMP-3 can also activate some pro-MMPs (e.g.

It is extremely important for pro-MMP-1 to produce MMP-1). MMP-11 is also called stromelysin-3, but it is often classified to other types because of its different structural sequence and substrate specificity.

2.4 Matrilysins

The structure of matrilysins lacks hemopexin domain. Matrilysin-1(MMP-7) and matrilysin-2(MMP-26) are included. In addition to degrading ECM, MMP-7 can process cell surface factors, such as pro TNF- α and E-cadherin. The matrilysin-2 (MMP-26) can also digest many ECM components.

2.5 Membrane-type MMPs(MT-MMPs)

There are six types of MT-MMPs, four types of MT-MMPs are type I transmembrane protein(MMP-14,-15, -16,-24), the other two types are GPI-ankyrin(MMP-17,-25). In addition to MT4-MMP, all other MT-MMPs can activate pro MMP-2. These enzymes can also digest some ECM components. MT1-MMP(MMP-14) has the activity to degrade collagen types I, II, III. MT5-MMP(MMP-24)

is mainly expressed in the cerebellum. MT6-MMP(MMP-25) is almost exclusively expressed in peripheral blood leucocytes, the original astrocytomas and the glioblastomas. However, it is not expressed in the meningiomas.

MT-MMPs have the following characteristics:

- a. Binding with cell membrane and provide the focus area for the decomposition and degradation of extracellular matrix protein;
- b. In the process of migration to the cell membrane, MT-MMPs are responsible for the cell activation of pro convertase pathway. Therefore, different from other types of MMPs, MT-MMPs have the proteolytic activity once it is inserted into the cell membrane.
- **c.** MT-MMPs have the substrate recognition sites of other types of MMPs and participate in an important pathway to activate other types of MMPs. It has been proved that MT1-MMP participates in the hydrolysis of MMP-2 and MMP-13.

2.6 Other types of MMPs

There are at least seven types of MMPs which are not mentioned in the prevous text. MMP-12 is mainly expressed in macrophages and may take part in the migration of macrophages. MMP-12 can digest elastin and some other proteins. MMP-19 is identified through the liver c-DNA cloning and exists as T-cell antigens in the rheumatoid arthritis patients. Enamelysin (MMP-20) is mainly expressed in newly formed enamel and can digest amelogenin. The function of MMP-22 is not clear. MMP-23 is mainly expressed in the regeneration tissue and Epilysin (MMP-28) is mainly expressed in keratinocytes^[7].

3. Influence factors to the activation of matrix metalloproteinases

The activation of MMPs is adjusted by three parts: gene transcription, zymogen activation and endogenous inhibitors.

3.1 Gene transcription

The genetic study has confirmed that there is promotor polymorphisms in MMP-1, -3, -9, -12 which affects the expression of MMPs gene^[8]. The expression of MMPs mRNA is affected by

a variety of chemical factors, such as ECM components, carcinogens, oncogene, neural hormones, cytokines and corticosteroids. For example, TNF- α affects the expression of MMPs gene by affecting the multi-cell system transcription factors (the latter is combined with the specific response elements in the MMPs gene enhancer);Tissue-type plasminogen activator (tPA) and urokinase can induce the expression of MMP-9^[9,10]. The signal transduction mechanism of MMPs' activation remains to be elucidated. Previous studies suggest that MAPK pathway may be relevant to MMPs gene expression, transcription factor (AP-1 and NF-KB) is closely relevant to activation of MMPs ^[10]. The signal transduction mechanisms in different types of cells or for different types of MMPs may be different. Mengshol et al^[11] have showed that p38, JNK and NF-KB are essential for IL-1 to induce

cartilage cells to produce MMP-13. For the production of MMP-1, p38 is still essential while JNK and NF-KB are not essential. Studies have also shown that some factors such as glucocorticoid and TGF- β , can inhibit MMP gene expression at the genetic level.

3.2 Zymogen activation

After translation and modification, the vast majority of MMPs mRNA is secreted to ECM in the form of zymogen. And it will be activated after hydrolysis of the propeptide domain. The activation mechanism found now includes stepwise activation, activation by MT-MMP, and cell activation. The initial activation of MMP is often associated with plasmin parenzyme, elastase, kallikrein. Among them, plasmin is considered to be the most powerful physiological activator in human body. In addition, SH-reagent (iodoacetic acid, HOCL, oxidized glutathione), denaturant (urea, SDS, NaSCN) and heat treatment have the ability to hydrolyze propeptide domain.

3.3 Endogenous inhibitors

The proteolytic activity of MMP is inhibited by non-specific and specific inhibitors. Nonspecific inhibitors include α 2-macroglobulin, α 1-antiprotease and BB-94(batimistat). Specific inhibitors are the tissue inhibitors of metalloproteinases(TIMPs). TIMPs, which are the coding proteins of multi-gene family, are the natural inhibitors of MMPs. The expression of TIMPs is regulated during development and tissue remodeling. To date, a total of four types of TIMPs have been found in vertebrates. They form the high affinity complexes with activated MMPs at the molar ratio of 1:1 and inhibit the degradation of ECM by blocking the catalysis domain of MMPs. TIMP-1 inhibits the activity of most MMPs, except for MT1-MMP and MMP-2; TIMP-2 inhibits majority of MMPs except MMP-9. In addition, TIMP-2 can form the complexes with MT1-MMP in the cell membrane which may have the regulation to activate the proteolytic activity of MMP-2; TIMP-3 inhibits MMP-1, -2, -3, -9, -13; TIMP-4 inhibits MMP-1, -3, -7, -9 and is highly expressed in human heart. Corresponding to MMPs, TIMPs play a negative role in the regulation of the ECM metabolism. They can prevent the activation of MMPs, inhibit the function and affect the extent of protein breakdown and duration of injury^[12]. In some excessive matrix degradation diseases, the imbalance between MMPs and TIMPs leads to a net increase in overall activity of MMPs^[13]. Although TIMPs play an important role in preventing excessive matrix degradation caused by MMPs, recent research shows that TIMP-1 and TIMP-2 are the multifunctional proteins with different biological functions. It has been reported that TIMP-1 and TIMP-2 demonstrate growth factor-like activity and inhibit the angiogenesis; while TIMP-3 is

148

www.intechopen.com

associated with apoptosis^[7].

TIMPs inhibit the MMPs activity by two steps. In the stage of zymogen activation, TIMP-2 can form a stable complex with pro MMP-2, as well as TIMP-1 and pro MMP-9. Therefore, they can impede the zymogen self-activation of pro MMP-1. In the stage of activated MMPs, both TIMP-1 and TIMP-2 can directly form a tight 1:1 complex with the activated MMPs and inhibit their activity. Naturally, the TIMP-MMP complexes may also be activated under certain conditions.

4. Role of matrix metalloproteinase in cerebral ischemia

In recent years, there is growing research interest on MMPs in central nervous system. In the normal central nervous system, MMP-2 and MMP-9 have been found in perivascular cells, brain vascular endothelial cells, astrocytes and microglias. The microglias in cultured rat can secrete MMP-9 once they are activated. Human hippocampal pyramidal cells can also synthesize MMP-9. It is reported that MMP-9 mRNA can be expressed in the developing mouse embryo brain, suggesting that MMP-9 is related to neurodevelopment. Studies also show that MMP-2 is related to the regeneration of axons. A growing number of scholars believe that MMPs plays an important role in the pathologic processes of central nervous system diseases^[5].

4.1 MMPs and pathogenesis of stroke

Atherosclerosis is one of the underlying vascular risk factors for developing cerebrovascular disease. It is reported that the gelatinase of MMPs plays a key role in the process of intimal injury and the formation of atherosclerotic lesions. The endothelial cells covered on the plaque have the activity to express MMP-2 and MMP-9. The over-expression of MMPs can dissolve collagens and significantly change the proportion of plaque composition which leads to the relative increase in lipid content and increases the plaque instability. Then the plaque cap thinningzs and splits, eventually leading to the occurrence of cerebral ischemia disease^[14,15]. Studies have shown that MMP-8 and MMP-1, MMP-12 might play a decisive role in maintaining the stability of atherosclerotic plaques^[16]. The rupture of atherosclerotic plaque partially depends on the activity and content of MMP-9^[17]. Another study reported that the focal increase of MMP-9 activity was an early warning of acute plaque rupture.

4.2 MMPs and ischemic brain injury

In the normal state, cerebral vascular endothelial cells can not express or only express a small amount of MMP-9^[18]. By the animal experiments, Fukuda^[19] confirmed that ischemic brain tissue could produce some active proteases and these proteases led to the rapid and significant degradation of microvessels. Another studies suggested that cerebral ischemia and reperfusion could induce the expression of MMPs. Especially, the activity of MMP-2 and MMP-9 would increase, which was related to cerebral microvascular permeability, blood-brain barrier (BBB) permeability, BBB damage, inflammatory cell invasion and cerebral edema^[6,12,13,20]. Gidday et al^[21] showed that MMP-9 in ischemic brain tissue played a pro-inflammatory role which helped neutrophil leukocytes migrate from the blood circulation into the tissue. At first, MMP-9 caused BBB damage, then MMP-9 contributed to the proteolysis of microvascular basement membrane and eventually led to neurological

damage. The possible mechanism may be the following: In the ischemic brain tissue, MMP-9 is mainly expressed in vascular endothelial cells, the increased MMP-9 may act on the tight junctions and basement membrane among the BBB endothelial cells which leads to the BBB damage, increase of permeability and vasogenic brain edema. The worst is the occurrence of herniation^[22]. In addition, the degradation of vascular basement membrane makes neutrophil leukocytes exudated to the brain tissue. MMP-2 and MMP-9 expressed by macrophages may contribute to their entry into the ischemic lesions and promote the wound healing after focal stroke ^[13]. However, MMP-9 expressed by neutrophil leukocytes may also contribute to the ischemic damage of brain. Zalewska et al^[23] proposed that activated MMP-9 may act on a certain link of the cell apoptosis cascade in the hippocampus CAI area after transient ischemic brain. Previous experimental studies also showed that MMP-9 could degrade the myelin basic protein of brain white matter and lead to the damage after ischemic brain^[24]. Early studies have showed that injecting MMP-2 to the rat brain can contribute to the opening of BBB. MMP-2 and MMP-9 destroy the capillary tight junctions and basement membrane by protein hydrolysis after brain ischemia, thus leading to vasogenic brain edema. In 1996, Rosenberg et al^[25] monitored of the expression of MMP-9 in rats with first onset cerebral infarction. They found the upregulated expression of MMP-9 4 hours after the occlusion of middle cerebral artery,. Within 12 hours and 24 hours, the expression of MMP-9 in the infarction site significantly increased which was consistent with the peak of the vasogenic brain edema, suggesting that MMP-9 played an important role in secondary brain damage and vasogenic edema. Montaner' et al^[26] also showed that the activity of MMP-9 abnormally increased in the early stage of stroke and pro-inflammatory response, while the activity of MMP-2 appeared in the repair phase of vascular regeneration. Autopsy results showed that one week after infarction, MMP-9 was expressed in neutrophils, and one week later the macrophages expressing matrilysin and MMP-2 were observed. Montaner et al also found that the level of plasma MMP-2 was higher in the patients with previous history of stroke. The above results suggested that MMP-9 plays an important role in secondary brain damage and vasogenic brain edema, while MMP-2 is involved in tissue repair and nerve regeneration.

Previous studies also suggested that MMP-9 may also related to the hemorrhage translation after tPA or urokinase thrombolysis^[24].

4.3 MMPs and cerebral ischemia reperfusion

Animal model studies showed that 3 hours and 48 hours of reperfusion occurred after cerebral ischemia, BBB was opened and the opening reached to the peak in the 48th hour; The first opening was related to the increased level of MMP-2 while the second opening occurred in the stage when the level of MMP-9 was significantly higher; the content of MMP-2 reached to the peak 5 days after reperfusion and the repair process began at the same time; the content of TIMP-1 significantly increased in the 48th hour while that of TIMP-2 increased to the maximum on day 5^[12]. These results suggested that reperfusion may affect MMPs and TIMPs, while MMPs and TIMPs promoted the reperfusion injury by complex ways. When the synthetic inhibitors of MMPs (BB-1101) was applicated to inhibit MMPs, both the BBB's first opening and cerebral edema after reperfusion were prevented, suggesting that BBB's opening and brain edema after reperfusion were related to MMPs.

150

5. Effect of inhibiting matrix metalloproteinase on cerebral ischemia

Matrix metalloproteinase inhibitors have been used to treat cancer metastasis. Currently, many MMP inhibitors are also used in experimental models of neurological diseases^[6], such as bacterial meningitis, cerebral infarction and experimental allergic meningitis. Studies have shown that the content of MMPs increased after cerebral ischemia and reperfusion, exogenous inhibitors could reduce the ischemic and reperfusion injury^[12,13]. Thus, MMPs may become a new potential target for stroke therapy and matrix metalloproteinase inhibitors can be used for treatment of cerebrovascular diseases^[2].

5.1 Matrix metalloproteinase inhibitors and cerebral ischemia

Previous studies have showed that the activated leukocytes could increase the reperfusion injury in the central nervous system. So, the drugs which can inhibit the leukocyte adhesion (including the intracellular adherence factor antibodies) have a neuroprotective effect. Matrix metalloproteinase inhibitors can combine with the divalent cation in vitro, inhibit the leukocytes function, and reduce the reperfusion injury.

Lee' et al^[3] have shown that MMP inhibitors could inhibit the MMP-9 production in brain after stimulated and parenchymal angiogenesis.Horstmann^[27] found that MMP-1 not only had the direct proteolytic capacity, but also played a role in the activation cascade of MMPs. It could crack collagen types I, II, III and be involved in the activation of MMP-2, MMP-9. Non-specific matrix metalloproteinase inhibitors are clustered in atherosclerosis tissues and inhibit the activity of MMP-9 in the carotid artery plaques, thereby stabilizing the easily broken atherosclerotic plaques. Matrix metalloproteinase inhibitors can also reduce the incidence of acute plaque rupture by reducing MMP-9 activity. The content of TIMPs is significantly higher in cerebral ischemia-reperfusion. They combine with the corresponding MMPs and prevent the activation of MMPs. They inhibit the function of MMPs, stabilize the ECM, significantly reduce the blood-brain barrier damage and the brain edema after ischemia^[12].

5.2 Matrix metalloproteinase inhibitors and therapies for cerebral ischemia

Previous studies suggested that MMP-9 monoclonal antibody may significantly reduce the infarction volume in a rat model of ischemia^[13]. Clinical studies have also confirmed that MMP-9 is related to the total infarction volume. It is reported that matrix metalloproteinase inhibitors and MMPs neutralizing antibodies can reduce the vasogenic brain edema and infarction volume^[26]. In addition, MMPs inhibitors are also effective in preventing atherosclerosis and the ischemic brain damage.

rt-PA and urokinase are effective drugs for acute ischemic stroke. However, the clinical application of thrombolytics is limited by the narrow therapeutic time window and the complication of bleeding after thrombolysis. But if we combine rt-PA or urokinase with the MMP inhibitors (BB-94 or doxycycline)in thrombolysis treatment, the incidence of hemorrhage and the amount of hemorrhage after thrombolysis may be decreased and the thrombolytic time window will be prolonged^[28,29]. This is because that if MMPs inhibitors are used before the application of rt-PA or urokinase, MMP-2, MMP-3 and MMP-9 (which have the potential damage for BBB) would be inhibited. BBB would be closed and the integrity of blood vessels would be maintained, thus might increase the safety of the thrombolytic therapy^[6,30].

6. Conclusion

MMPs may play an important role in the brain ischemia and reperfusion by degrading the ECM and destroying the blood -brain barrier which can lead to the vasogenic brain edema and secondary brain injury. The matrix metalloproteinase inhibitors and MMPs neutralizing antibodies can reduce the vasogenic brain edema and infarction volume. In addition, MMPs inhibitors are also effective in preventing atherosclerosis and the ischemic brain damage. Thus, MMPs may become a new potential target for stroke therapy and matrix metalloproteinase inhibitors can be used for the treatment of cerebrovascular diseases.

7. References


- [1] Sheppeck JE 2nd, Tebben A, Gilmore JL, et al.A molecular modeling analysis of novel non-hydroxamate inhibitors of TAGE.Bioorg Med Chem Leit,2007,17:1408-1412.
- [2] Hashimoto T, Wen G, Lawton MT, et al. Abnormal expression of matrix metalloproteinases and tissue inhibitors of metalloproteinases in brain arteriovenous malformations. Stroke,2003,34:925-931.
- [3] Lee CZ, Xu B, Hashimoto T, er al. Doxycycline suppresses cerebral matrix metalloproteinase-9 and angiogenesis induced by focal hyperstimulation of vascular endothelial growth factor in a mouse model.Stroke,2004,35:1715-1719.
- [4] Newby AC.Dual Role of Matrix Metalloproteinases(Matrixins) in Intimal Thickening and Atherosclerotic Plaque Rupture.Physiol Rev,2005,85:1-31.
- [5] Axiso B, Loftus I, Ross Naylor A. Prospective, randomized, double-blind trial investigating the effect of Doxycycline on matrix metalloproteinase expression within delayed carotid plaques(J).Stroke,2002,33:2858-2865.
- [6] Pfefferkorn T, Rosenberg GA. Closure of the blood-brain barrier by matrix metalloproteinase inhibition reduces rtPA-mediated mortality in cerebral ischemia with delayed reperfusion.Stroke,2003,34:2025-2030.
- [7] Visse R, Nagase H. Matrix metalloproteinases and tissue inhibitors of metalloproteinases: structure, function and biochemistry.Cire Res,2003,92:827-839.
- [8] Ye S. Polymorphism in matrix metalloproteinase gene promoters:implication in regulation of gene expression and susceptibility of various diseases.Matrix Biol,2000,19:623-629.
- [9] LI Jie,LI Guang-qin,JIANG Guo-hui.Effects of rt-PA and urokinase on expression of MMP-9 and blood-brain barrier permeability in thrombolytic therapy.Chin J Contemp Neurol Neurosurg,2009,9:445-449.
- [10] ZHANG Shu-tang,LI Guangqin.The effect of urokinase on the expression of MMP-9 in human endothelial cell and its mechanism.Journal of Chongqing Medical University,2009,34:269-272.
- [11] Mengshol JA, Vincenti MP, Coon CI, et al. Interleukin-1 induction of collagenase3 (matrix metelloproteinase 13) gene expression in chondrocytes requires p38,c-Jun N-terminal kinase,and nuclear factor kappaB: differential regulation of collagenase 1 and collagenase 3.Arthritis Rheum,2000,43:801-811.
- [12] Rosenberg GA, Estrada EY, Dencoff JE. Matrix metalloproteinases and TIMPs are associated with blood-brain barrier opening after reperfusion in rat brain. Stroke,1998,29:2189-2195.

- [13] Romanic AM, White RF, Arleth AJ, er al. Matrix metalloproteinases expression increases after cerebral focal ischemia in rats: inhibition of matrix metalloproteinase-9 reduces infarct size. Stroke,1998,29:1020-1030.
- [14] Morgan AR, Rerkasem K, Gallagher PJ, et al. Differences in matrix metalloproteinase-1 and matrix metalloproteinase-12 transcript levels among carotid atherosclerotic plaques with different histopathological characteristics. Stroke,2004,35:1310-1315.
- [15] Fatar M, Stroick M, Griebe M, er al. Matrix metalloproteinases in cerebrovascular diseases. Cerebrovasc Dis,2005,20:141-148.
- [16] Molloy KJ, Thompson MM, Jones JL, et al. Unstable carotid plaques exhibit raised matrix metalloproteinase-8 activity. Circulation,2004,110:337-343.
- [17] Heo JH, Kim SH, Lee KY, et al. Increase in plasma matrix metalloproteinase-9 in acute stroke patients with thrombolysis failure. Stroke,2003,34:e48-50.
- [18] Guzener ML, Gveric D, Strand C, et al. The expression of tissue-type plasminogen activator, matrix metalloproteinases and endogenous inhibitors in the central nervous system in multiple sclerosis: comparision of stages in lesion evulution.J Neuopathol Exp Neurol, 1996, 55:1194.
- [19] Fukuda S, Fini CA, Mabuchi T, et al. Focal cerebral ischemia induces active proteases that degrade microvascular matrix(J). Stroke,2004,35(4):998-1004.
- [20] Gasche Y, Fujimura M, Morita-Fujimura Y, et al. Early appearance of activated matrix metalloproteinase-9 after focal cerebral ischemia in mice: a possible role in bloodbrain barrier dysfunction. J Cereb Blood Flow Metab,1999,19:1020-1028.
- [21] Giddy JM, Gasche YG, Copin JC, et al. Leukocyte derived matrix metalloproteinase-9 mediates blood-brain barrier breakdown and is proinflammatory following transient focal cerebral ischemia. Am J Physiol Heart Circ Physiol,2005,289:H558-H568.
- [22] Eldrup N, Gronholdt ML, Sillesen H, et al. Elevated matrix metalloproteinase-9 associated with stroke or cardiovascular death in patients with carotid stenosis. Circulation,2006,114:1847-1854.
- [23] Zalewska T, Ziemka-Nalecz M, Sarnowska A, et al. Involvement of MMPs in delayed neuronal death after global ischmia. Acta Neurobiol Expe,2002,62:53-61.
- [24] LI Jie,LI Guang-qin.Effects of urokinase thrombolysis on matrix metalloproteinase-9, blood-brain barrier permeability and intracerebral hemorrhage in a rat modal of middle cerebral artery thromboembolism.Int J Cerebrovasc Dis,2009,17:668-672.
- [25] Rosenberg GA, Navratil M, Borane F, et al. Proteolytic cascade enzymes increase in focal cerebral ischemia in rat. J Cereb Blood Flow Metab,1996,16:360-366.
- [26] Montaner J, Alvarez-Sabin J, Molina C, et al. Matrix metalloproteinase expression after human cardioembolic stroke: temporal profile and relation to neurological impairment. Stroke,2001,32:1759-1766.
- [27] Horstmann S, Kalb P, Koziol J, et al. Profiles of matrix metalloproteinases: their inhibitors, and Jaminin in stroke patients: infulence of different therapies. Stroke,2003,34:2165-2170.
- [28] Lapchak PA,Chapman DF,Zivin JA. Metalloproteinase inhibition reduces thrombolytic (tissue plasminogen activator)-induced hemorrhage after thromboembolic stroke.Stroke,2000,31:3034-3040

- [29] JIANG Guo-hui,LI Guang-qin,LI Jie.Inhibtion of matrix metalloproteinase-9 can prolonging urokinase thrombolytic time window after embolic stroke in rats.Journal of Brain and Nervous Diseases,2009,17:354-358
- [30] Montaner J, Molina CA, Monasterio J, et al. Matrix metalloproteinase-9 pretreatment level predicts intracranial hemorrhagic complications after thrombolysis in human stroke. Circulation,2003,107:598-603.







Advances in the Preclinical Study of Ischemic Stroke Edited by Dr. Maurizio Balestrino

ISBN 978-953-51-0290-8 Hard cover, 530 pages Publisher InTech Published online 16, March, 2012 Published in print edition March, 2012

This book reports innovations in the preclinical study of stroke, including - novel tools and findings in animal models of stroke, - novel biochemical mechanisms through which ischemic damage may be both generated and limited, - novel pathways to neuroprotection. Although hypothermia has been so far the sole "neuroprotection" treatment that has survived the translation from preclinical to clinical studies, progress in both preclinical studies and in the design of clinical trials will hopefully provide more and better treatments for ischemic stroke. This book aims at providing the preclinical scientist with innovative knowledge and tools to investigate novel mechanisms of, and treatments for, ischemic brain damage.

How to reference

In order to correctly reference this scholarly work, feel free to copy and paste the following:

Wan Yang and Guangqin Li (2012). The Matrix Metalloproteinases and Cerebral Ischemia, Advances in the Preclinical Study of Ischemic Stroke, Dr. Maurizio Balestrino (Ed.), ISBN: 978-953-51-0290-8, InTech, Available from: http://www.intechopen.com/books/advances-in-the-preclinical-study-of-ischemic-stroke/matrix-metalloproteinase-and-cerebral-ischemia

INTECH

open science | open minds

InTech Europe

University Campus STeP Ri Slavka Krautzeka 83/A 51000 Rijeka, Croatia Phone: +385 (51) 770 447 Fax: +385 (51) 686 166 www.intechopen.com

InTech China

Unit 405, Office Block, Hotel Equatorial Shanghai No.65, Yan An Road (West), Shanghai, 200040, China 中国上海市延安西路65号上海国际贵都大饭店办公楼405单元 Phone: +86-21-62489820 Fax: +86-21-62489821 © 2012 The Author(s). Licensee IntechOpen. This is an open access article distributed under the terms of the <u>Creative Commons Attribution 3.0</u> <u>License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original work is properly cited.

IntechOpen

IntechOpen